

Heater BB.WARMER



BB.WARMER with heating trough uses heating rods used to support the treatment of pediatric patients during emergency resuscitation or in other medical care environments that require a different heating temperature. The outstanding technical features of BB.WARMER include:

- Automatically adjust the heating energy based on the baby's body temperature change and the ambient temperature
- Automatic warning mode brings absolute safety in treatment.
- Simple operation and easy to use even for first time operators.



Specifications

Measure	<ul style="list-style-type: none"> - Display resolution: 0.1⁰C - Measurement resolution: 0.125⁰C (1/8⁰C) - Measurement parameters: ambient temperature, probe (sensor) temperature
Measuring range	<ul style="list-style-type: none"> - Measurement temperature of Heating:-55 ÷ 125⁰ C - Operating temperature setting: 32÷ 39⁰C
Control	- Control the heating power automatically proportionally based on the set temperature and the ambient temperature, ensuring safety by the probe temperature.
Size	<ul style="list-style-type: none"> - Heating tray: 460 x 222 x 112 mm (LxWxH) - Heating foot dimensions: 610 x 480 mm (LxW)
Height adjustment	- From 1600 to 1800 mm (Adjustable range 200 mm)
Heating rod	- 220VAC – 400 W
Heating rod duration	- Minimum 2000 hours of use
Heating area	800 x 600 mm (LxW) at 700 mm distance from heater to pediatric patient
Wheel	4 swivel castors, diameter 65 mm
Adjust the projection angle of the heater	± 40°
Weight	16 kg
Infusion hanger	150 mm
Work environment	<ul style="list-style-type: none"> - Temperature: From 0°C to 40⁰C - Humidity: from 0% to 98%, no condensation
Transport and storage requirements	<ul style="list-style-type: none"> - Temperature: from 0°C -50°C - Humidity: from 0% to 99%, no condensation
Noise	≤ 20 dB (or 15 dBA - level corrected for the same frequency to match human ear characteristics)
Power Consumption	≤ 500W
Voltage	185 to 240VAC, 47-63Hz
Optional accessories	Temperature sensor

Products by Manufacturer: VINH DAT .LLC - BB TIME

Room 810, 8th floor, JSC34 Building, Alley 164 Khat Duy Tien, Thanh Xuan, Ha Noi

Tel: 04-222 50797 E-

mail : sale@bbtime.vn hoặc vn.bbtime@gmail.com

Fax : 04-222 50643

Website : www.bbtime.vn